



## LAUREL BRIDGE


### NEWS RELEASE

Press Contact: Susan Reagan  
Tel: 302-453-0222  
E-mail: [susan@laurelbridge.com](mailto:susan@laurelbridge.com)

### Laurel Bridge Software Announces a DICOM Archive Migration Toolkit

**Newark, DE (May 1, 2005)** – Laurel Bridge Software announced the release of the Laurel Bridge Migration Toolkit (LBMT). The LBMT is a collection of applications that provides support for companies transitioning from one DICOM compatible image archive to another. The centerpiece of the LBMT is the Laurel Bridge Migration Controller (LBMC) which manages all the migration related activities.

The LBMC provides automated facilities and services that allow migration of studies, series and SOP instance data from one DICOM 3.0 compatible physical archive to another. The Migration Controller application supervises and schedules the transfer, plus provides verification to ensure that a migration is complete and accurate. Migration history is logged and performance statistics are gathered; various levels of reports are provided to summarize the status and progress of a migration. Additional cost options are available to provide site-specific mapping and correction of DICOM elements that may be required when moving to a new archive. A typical migration is illustrated in the diagram below.



**Laurel Bridge Software, Inc.** is a Delaware-based corporation, specializing in providing software tools and utilities for use by software professionals, especially in areas related to networking and controlling medical devices. Engineers from Laurel Bridge and its sister company, Blair Computing Systems, Inc., have been developing medical imaging systems for more than a decade. Together they have a wealth of experience developing PACS networks, modalities, workstations/viewers, archives and protocol or image converter boxes of all shapes and sizes, working with companies like AGFA, Codonics, Direct Radiography Corporation, DuPont, Hologic, LORAD, and others. Contact: Laurel Bridge Software, Inc., 160 East Main Street, Newark, DE 19711 USA, Tel: 302-453-0222, Fax: 302-453-9480, E-mail: [info@laurelbridge.com](mailto:info@laurelbridge.com), or visit [www.laurelbridge.com](http://www.laurelbridge.com).

DICOM® is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communications of medical information. Laurel Bridge, its logo, DPA, and DCF are trademarks of Laurel Bridge Software, Inc.