



FOR IMMEDIATE RELEASE

Contact: Greg Muller
(302) 453-0222

Greg.Muller@LaurelBridge.com

Laurel Bridge Software Announces ‘Immediate Mode’ Processing for Compass DICOM Router

Newark, DE (February 4, 2013) – Laurel Bridge announces release 2.1.35 of Compass™, its sophisticated, rules-based DICOM® Store and Forward Routing application for DICOM store jobs. Selected feature improvements include:

- Widely requested, the Immediate Mode processing option allows Compass to promptly route and forward images while an incoming study is still being received. By reducing study transfer latency, this feature significantly improves image workflow efficiency.
- TLS/SSL secure communication enhancements for use with Immediate Mode transfers.
- Improved logging for destination configuration and for Automatic Heartbeat Sensing network status monitoring.
- Expanded rules for selectively matching text-based DICOM tags used in tag-morphing.
- Hot folder operational improvements including persisting file ingest failures to allow operator intervention.

Laurel Bridge continues to aggressively solicit and incorporate customer feedback and recommendations into the Compass product.

About Laurel Bridge Software

Laurel Bridge Software provides solutions, tools, and services to simplify medical image workflows and DICOM networking. Application solutions route, store, filter, and migrate medical image data securely and reliably. DICOM software development toolkits (SDKs) provide application developers the ability to integrate DICOM functionality into their products. Laurel Bridge also offers a full complement of integration, configuration, application development, and technical support services.

For more information, contact Greg Muller at info@LaurelBridge.com, or visit www.LaurelBridge.com.

DICOM® is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to Digital Imaging Communication of Medicine. Laurel Bridge, its logo, Simplifying DICOM, and its product names, including DCF, Exodus, Switchboard, Compass and Power Tools, are trademarks of Laurel Bridge Software, Inc.